







How We Compare

| |  | Discontinued  | Discontinued  |  |
|-----------------------|---|---|--|---|
| Prescription Needed | No | No | Yes | Yes |
| Adjustable Pressure | Yes | No | No | No |
| Expiratory Resistance | Variable (Provides a wide range) | Fixed (Equivalent to Optipillows with side vent 60% closed) | Fixed (Equivalent to Optipillows with side vent 75% closed) | Fixed (Equivalent to Optipillows with side vent 60% closed) |
| Reusable | Yes | No | No | Yes |
| Application | Soft nasal pillows with headgear | Adhesives | Adhesives | Nasal Insert (headgear optional) |

Optipillows EPAP Mask is the only EPAP device that allows users to adjust the pressure by opening or closing the valve (range is between 1.5 cmH₂O and 15 cmH₂O). Optipillows is reusable and cost effective. The soft nasal pillows are comfortable and rest gently against the face to create a perfect seal. Optipillows is lightweight and perfect for travel.

| Peak Expiratory Pressure generated by the EPAP Mask | | |
|--|------------------------------------|-------------------------------------|
| With Side Vent Fully Open | With Side Vent 50% Closed | With Side Vent 75% Closed |
| 1 to 3 cmH ₂ O pressure | 3 to 7 cmH ₂ O pressure | 8 to 14 cmH ₂ O pressure |

Email: michelle.optipillows@gmail.com Phone: Michelle 360-921-6480
www.Optipillows.com